



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 8141**



## **WORK SURFACE, LABCONCO (3 ft, Stainless Steel)**


**W 36" × D 29" × H 0.8"**

**DESCRIPTION:** This stainless steel work surface helps to contain spills that may happen while using the enclosure.

### **TECHNICAL SPECIFICATIONS:**

Dimension, W × D × H (cm / in)	91.4 × 73.7 × 2.0 / 36 × 29 × 0.8
Shipping Weight (Kg / lb)	45.4 / 100
Material	Stainless Steel

### **ADDITIONAL INFORMATION:**

-  The enclosure should be sealed at the work surface to prevent spilled materials from collecting under the walls