

8178-01



Flip-Off Seal, Plastic w/Aluminum Base, Dark Grey, 20 mm, pk/1000 (WHEATON)

A Flip-Off seal is your drug product's first line of protection. It helps to prevent contamination from getting into the vial and seals the parenteral package. Further, the seal helps to protect the injection site on the vial stopper. Recommended seals consist of an aluminum shell and a plastic Flip-Off button.

The above information is a direct transcription of information provided to MEDISCA® from the Specification Sheet and/or Certificate of Conformity provided by the manufacturer/supplier.