

PRODUCT NO.: 8209



## Bottle w/Tip, 2oz, Translucent, LDPE, 20-410

**Description:** These 2oz Enema Bottles are manufactured from an especially pliable, low-density polyethylene. They provide quick and total dispensing of preparations from bottle to patient. The bottles are translucent so you can immediately see the fluid level contained within them. The distinctly lubricated enema bottle tips allow for easy insertion and removal from the patient.

### Features:

- Pliable Low Density Polyethylene
- Natural Color for Easy Fluid Identification
- Enema Tips are pre-lubricated

**\*Note:** Not intended as medical device

USA | Canada | Australia | International

medisca.com