



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8224



BEAKER, GRADUATED GLASS (2000 mL, Heavy Duty)

Graduation Range: 200–1800 mL, Graduation Interval: 200 mL, Dual Scale

DESCRIPTION: This heavy duty beaker with spout is specifically designed to offer good mechanical strength under harsh conditions, such as mechanized washing operations. For convenience, it has a double graduated metric scale to indicate approximate content and has an extra-large marking spot. It also conforms to ASTM E-960 standards.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	2000
Material	Type I Class A - Low Expansion Borosilicate glass
Graduation Range (mL)	200 – 1800
Graduation Interval (mL)	200
Accuracy	±5% of full capacity (±100 mL)
Approximate O.D. × H (mm)	131 × 180