



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 8454



BOSTON ROUND GLASS BOTTLE w/DROPPER (Amber, 1 oz / 30 mL, 20-400)
Type III Glass, Glass Dropper w/Black Cap

DESCRIPTION: ASTM E438, Type III amber glass protects light-sensitive chemicals. It is supplied with 20-400 GPI screw cap, bulb and dropper assembled on the bottle. The dropper is USP Type III glass. There are no graduation markings on the dropper.