



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 8459**



**PLASTIC BOTTLE w/FLIP-TOP CAP (Natural, 4 oz / 120 mL, 24 – 410)**  
HDPE Plastic, White PE Unlined Cap

**DESCRIPTION:** This Polyethylene (PE) Plastic Flip Top Dispensing Bottle is ideal for dispensing small amounts of liquid.

**TECHNICAL SPECIFICATIONS:**

Capacity / Volume (mL / oz)	120 / 4
Neck Finish	24 – 410
Size, Ø x H (mm)	47 x 124
Material	HDPE (High Density Polyethylene)
Color	Natural
Bottle Style	Cylinder
Attached Cap	White PE (Polyethylene) Unlined Flip Top Cap