



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO. 8670






4 ft TELESCOPING BASE STAND WITH FIXED FEET (LABCONCO)


For XPert / RXPert Filtered Balance Systems

DESCRIPTION: Telescoping Base Stands with Fixed Feet may be adjusted to eight height positions prior to installation. Telescoping legs are adjustable in 1" increments from 27.5" to 34.5". Durable 1.75" tubular, epoxy-coated steel frame. Rear cross bar supports on all stands except 2' width.

FEATURES:

-  Weight: 110.0 lbs
-  Weight metric: 50.0 kg
-  Dimensions: 48.0" w x 29.0" d
-  Dimensions metric: 121.9 x 73.7 cm
-  Conformance: ADA
-  Maximum Stand Height: 34.5" (85.1 cm)
-  Minimum Stand Height: 27.5" (69.9 cm)
-  Nominal Width: 4'

RELEVANT PRODUCTS:

-  Applicable Products: Basic Hoods, Paramount Ductless Enclosures, Protector Hoods, Protector Work Stations, Protector XVS Ventilation Stations, Purifier Class I Enclosures, Purifier Filtered PCR Enclosures, Purifier HEPA Filtered Enclosures, Purifier Horizontal Clean Benches, Purifier Vertical Clean Benches, XPert Balance Enclosures, XPert Filtered Balance Stations, Xpert & RXPert Filtered Balance Systems, XPert Nano Enclosures