



PRODUCT NO.: 8759



NEUTRALIZING ACID RINSE (Labconco)

34oz, Reduces Alkaline Detergent Residue, 6 Uses, For Labconco SteamScrubbers

DESCRIPTION: Specially formulated for rinse aid reservoir in Labconco glassware washers. Mild acid corrects pH reducing alkaline detergent residue. Removes hard water deposits, mineral deposits, and accumulated lime scale that can cloud glassware. Rinses completely. Fills reservoir approximately 6 times. 34 ounces. LabSolutions Powder Detergent, Low-Foaming Liquid Detergent, and Neutralizing Acid Rinse are specially formulated for Labconco SteamScrubbers and FlaskScrubbers. Each removes a wide range of laboratory contaminants and rinses completely.

TECHNICAL SPECIFICATIONS:

- Weight: 2.7 lbs
- Weight metric: 1.2 kg
- Number of Loads: Not applicable

ADDITIONAL INFORMATION:

Applicable Products:

- 7496-01: SteamScrubber, Freestanding
- 8654-01: SteamScrubber, Undercounter