



PRODUCT NO.: 8763



UTENSIL BASKET (Labconco)

8"W X 8"D X 4"H, Stainless Steel, For Labconco SteamScrubbers

DESCRIPTION: For miscellaneous labware such as stoppers, spatulas and other small lightweight items that will not pass through No.4 mesh sides and bottom. Retainer top is included. Type 304 stainless steel. 8.8" w x 8.0" d x 4.0" h (20 x 20 x 10). Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

TECHNICAL SPECIFICATIONS:

- Weight: 1.5 lbs
- Weight metric: 0.7 kg
- Dimensions: 8"w x 8"d x 4"h
- Applicable Rack Types: Spindle, Standard
- Rack Location: Lower, Upper

ADDITIONAL INFORMATION:

Applicable Products:

- 7496-01: Freestanding SteamScrubber
- 8654-01: Undercounter SteamScrubber