PRODUCT NO.: 8794





Adapters, Press-In Bottle, for 28mm Bottle

DESCRIPTION: Press-in bottle adapters allow for increased dosing accuracy and fit securely in any bottle with the appropriate size opening. The Press-in Bottle Adapters allow a user to easily draw liquid with an oral dispenser, ensuring accurate dosing while avoiding spills. The adapters fit flush with the bottle opening so that the original cap can be placed on the bottle.

FEATURES:

- Clear color
- FDA approved plastic material
- Made up of LDPE

Which Press-In adapter will fit my bottle?			
Find your bottle by closure code, outer thread diameter or inner neck diameter below.			
Adapter Stock Code	Outer Thread Diameter	Inner Neck Diameter (ID) ⁽²⁾	<u>Height</u>
8796	16.26 mm	11.18 mm	14.78 mm
8795	20.42 mm	15.14 mm	15.41 mm
8794	22.62 mm	17.8 mm	15.03 mm
6349	28.05	22.66 mm	15.38 mm
Note: MEDISCA cannot guarantee full standardization of all bottles on the market. Please test your adapter for final confirmation of fit.			

⁽¹⁾ xxx refers to the thread finish, such as 400, 410 or 415.

⁽²⁾ ID values are approximate.

USA | CANADA | AUSTRALIA PACIFIC | EUROPE WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.