



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING




PRODUCT NO.: 8795



Adapters, Press-In Bottle, for 24mm Bottle

DESCRIPTION: Press-in bottle adapters allow for increased dosing accuracy and fit securely in any bottle with the appropriate size opening. The Press-in Bottle Adapters allow a user to easily draw liquid with an oral dispenser, ensuring accurate dosing while avoiding spills. The adapters fit flush with the bottle opening so that the original cap can be placed on the bottle.

FEATURES:

-  Clear color
-  FDA approved plastic material
-  Made up of LDPE

Which Press-In adapter will fit my bottle?

Find your bottle by closure code, outer thread diameter or inner neck diameter below.

Adapter Stock Code	Outer Thread Diameter	Inner Neck Diameter (ID) ⁽²⁾	Height
8796	16.26 mm	11.18 mm	14.78 mm
8795	20.42 mm	15.14 mm	15.41 mm
8794	22.62 mm	17.8 mm	15.03 mm
6349	28.05	22.66 mm	15.38 mm

Note: MEDISCA cannot guarantee full standardization of all bottles on the market. Please test your adapter for final confirmation of fit.

(1) xxx refers to the thread finish, such as 400, 410 or 415.

(2) ID values are approximate.

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM