

PRODUCT NO.: 8896



PLATE STAND (7.8"W X 12.4"D X 22"H)

Slip-Resistant Foil, Maximum Load: 5Kg, Supporting Rod: 16mm Diameter

DESCRIPTION: Plate stand to be used to secure T25 Homogenizer in place during operation.

FEATURES:

Diameter of support rod: 16 mm

Max. load: 5 kg

Dimensions (W x D x H): 7.8 x 12.4 x 22 in / 200 x 316 x 560 mm

ADDITIONAL INFORMATION:

Applicable Products:

9966-01: T25 Digital Homogenizer

8294-01: Boss Head Clamp

7393-01: Strap Clamp

Choose ONE Dispersing Element:

7970-01 Dispersing Element 1-50ml

7398-01 Dispersing Element 10-1500ml

8821-01 Dispersing Element 100-2000ml

8990-01 Dispersing Element 25-1500ml