

8901-01






Press-in Bottle Adapter, Plastic, Natural, 28 mm, pk/25 (E.P.S. Inc)

Note: Bottle and syringe not included.

Just press the adapter into the bottle's opening, attach your oral syringe, and withdraw your dose of medication from the bottle. It's that simple!

Features:

-  Bottle Can be Capped with Adapter Inserted
-  Manufactured From Low-Density Polyethylene (LDPE)
-  This size fits comfortably into most bottles with necks in its inner diameter range of 21 mm to 23 mm

The above information is a direct transcription of information provided to MEDISCA[®] from the Specification Sheet and/or Certificate of Conformity provided by the manufacturer/supplier. .