



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

## PRODUCT NO. 8963










### 4 ft TELESCOPING BASE STAND WITH CASTER (LABCONCO)


For XPert / RXPert Filtered Balance Systems

**DESCRIPTION:** Telescoping Base Stands with Casters may be adjusted to eight height positions prior to installation. Telescoping legs are adjustable in 1" increments from 27.5" to 34.5". Casters provide mobility for short distances, such as to gain access to the area behind the stand for service or maintenance. Toe-locking, non-marking polyurethane casters, 5" diameter, have bearings for quiet operation. Durable 1.75" tubular, epoxy-coated steel frame. Rear cross bar supports on all stands except 2' width.

#### FEATURES:

-  Weight: 84.0 lbs
-  Weight metric: 38.0 kg
-  Dimensions: 48.0" w x 40.0" d
-  Dimensions metric: 121.9 x 101.6 cm
-  Conformance: ADA
-  Maximum Stand Height: 33.5" (85.1 cm)
-  Minimum Stand Height: 27.5" (69.9 cm)
-  Nominal Width: 4'

#### RELEVANT PRODUCTS:

-  Applicable Products: Paramount Ductless Enclosures, Protector Work Stations, Protector XVS Ventilation Stations, Purifier Class I Enclosures, Purifier Filtered PCR Enclosures, Purifier HEPA Filtered Enclosures, Purifier Vertical Clean Benches, XPert Balance Enclosures, Xpert and RXPert Filtered Balance Systems, XPert Nano Enclosures