





8992-01



### **Octagonal Magnetic Stir Bar, White, Teflon w/Alnico-V Magnet, 2" x 5/16"**

Extra strength Teflon stirring bars molded with alnico-V magnets permanently coated with white PTFE. Inert to almost all chemical reagents. Can be used with temperatures up to 250°C. Excellent electrical properties, good thermal stability and very low static and dynamic coefficients of friction.

#### **Technical Specifications:**

-  Dimensions: 2" x 5/16"
-  Color: White
-  Material: Teflon coated with PTFE w/Alnico-V Magnet
-  Autoclavable

*The above information is a direct transcription of information provided to MEDISCA® from the Specification Sheet and/or Certificate of Conformity provided by the manufacturer/supplier.*