



BEAKER, GRADUATED GLASS (250 mL, Heavy Duty)

Graduation Range: 25–200 mL, Graduation Interval: 25 mL, Dual Scale

DESCRIPTION: This heavy duty beaker with spout is specifically designed to offer good mechanical strength under harsh conditions, such as mechanized washing operations. For convenience, this beaker has a double graduated metric scale to indicate approximate content and has an extra-large marking spot. It also conforms to ASTM E-960 standards.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	250
Material	Type I Class A - Low Expansion Borosilicate glass
Graduation Range (mL)	25 – 200
Graduation Interval (mL)	25
Accuracy	±5% of full capacity (±12.5 mL)
Approximate O.D. × H (mm)	68 × 91