



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9256-01



BOSTON SQUARE GLASS BOTTLE w/SCREW CAP
1000mL, 45mm, Amber Glass w/Blue Screw Cap

DESCRIPTION: Amber Boston Square Bottles protect contents from UV rays and are ideal for light sensitive products. These general use bottles are perfect for liquids. These bottles were designed for use with Socorex Calibrex™ 520 and 521 bottle-top dispensers.

FEATURES:

- Style: Square shaped amber glass reagent bottle
- Size, W x H (mm): 95 x 216
- Color: Amber, USP Type III
- Material: Glass
- Neck: 45mm
- Attached Cap: Blue Screw Cap
- Capacity/Volume: 1000mL

ADDITIONAL INFORMATION:

Applicable Products:

- 8961-01: Bottle-Top Dispenser, 1-10mL