



PRODUCT NO.: 9576



DROPTAINER W/TIP & CAP (3mL, STERILE)

2 per Pack, Plastic, Natural

DESCRIPTION: The sterile ophthalmic droptainer bottles are ideal when extemporaneous or small batches of sterile eye medications need to be prepared. Each droptainer blister contains 2 bottles, 2 tips and 2 caps. All droptainer bottles are manufactured from a special zinc stearate-free resin. The special plastic minimizes the risk of particulate forming in solution with certain ophthalmic preparations.