



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9615



HALF HEIGHT SIEVE (Ø 8")

#80 Mesh Size, Brass

DESCRIPTION: This sieve is durable and ideal for sieving applications where a medium sample size is required.

TECHNICAL SPECIFICATIONS:

Dimension, Ø × H (cm / in)		20.3 × 2.54 / 8 × 1
Mesh Size (US)		#80
Sieve Opening (µm / in)		180 / 0.0071
Material	Frame	Brass
	Wire Cloth	Stainless Steel
Conformance		ASTM E11, ISO 3310