



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9649



Spatula, Stainless Steel, w/Wood Handle, 100mm (4")

DESCRIPTION: This spatula is used to transfer solid ingredients and any other functions as required.

FEATURES:

- ✔ Multiple sizes allow for varied laboratory applications
- ✔ Made from ASTM stainless steel
- ✔ Non-sticking, corrosion-free standard blade
- ✔ Riveted, smooth hardwood handle designed and engineered to fit comfortably in your hand
- ✔ Electronically passivated for corrosion resistance
- ✔ Autoclavable
- ✔ Latex-free

TECHNICAL SPECIFICATIONS:

Blade Material	Stainless Steel
Handle Type	Riveted
Blade Length (cm/in)	10.16/4
Handle Material	Hardwood

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM