



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9662



**PLASTIC OINTMENT TUBE w/CAP (4 oz)**  
White, 4.75 mm Opening, #16 Neck Size

**DESCRIPTION:** This plastic ointment tube is lightweight yet durable to withstand squeezing and folding without developing cracks or leaks. The tube is made from Polyethylene extrusion grade resin blend.

**TECHNICAL SPECIFICATIONS:**

Dimensions w/Cap, Ø × L (cm / in)	3.49 × 17.15 / 1.375 × 6.750
Capacity (mL / oz)	118.3 / 4
Material	PE
Exterior color	White
Orifice (mm / in)	4.7 / 0.187