

PRODUCT NO.: 9758



pH-INDICATOR STRIP, NON BLEEDING

Measuring Range: 0 - 14

DESCRIPTION: This transparent pH-indicator strip for turbid solutions eliminates the need for sample preparation steps, such as filtration or clarification. The pH reaction zone is backed by a transparent carrier film, so the pH can be easily read on the reverse side of the strip.

The special indicator dyes are covalently bound to the reagent paper which effectively prevents the indicator from bleeding, and allows the strip to be left in the measurement medium indefinitely – without contaminating the sample.

FEATURES:

- 2 Non-bleeding strip prevents contamination of the medium
- 7 Transparent strip for clear results in lightly colored or turbid liquids
- Quick and easy method without sample preparation
- Brilliant color scales ensure reliable results
- SafetyEdge box for maximum security and convenience

ADDITIONAL INFORMATION:

- Store at 10 25 °C to maintain top condition for 3 5 years
- Protect from light and moisture
- Close box immediately after removing each strip or paper

Page 1 of 1

📰 1.800.932.1039 | 🖬 1.800.665.6334 | 🛅 1300 786 392 | WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.