



PRODUCT NO.: 9850



CYLINDER, GRADUATED GLASS TYPE 1 CLASS B (50 mL)

Graduated Range: 2.0 – 50.0 mL, Graduation Interval: 1.0 mL, Pharmaceutical Grade, Double Scale, w/Plastic Base & Protective Hexagonal Bumper Guard

DESCRIPTION: This cylinder is made of two parts: a strong, accurate glass graduate and a detachable, sturdy plastic hexagonal base which absorbs impacts and reduces breakage. The hexagonal base provides extra stability and all markings are durable, white enamel. Comes with a protective, polyethylene, hexagonal bumper guard.

This cylinder meets the ASTM E-1272 Class B, Type I specifications for classification and tolerances.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	50
Approximate O.D. x Height (mm)	24 x 212
Graduation Range (mL)	2.0 – 50.0
Graduation Interval (mL)	1.0
Tolerance (mL)	± 0.40