



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9853**



**CYLINDER, GRADUATED GLASS TYPE I, CLASS A (500 mL)**

Graduated Range: 25.0 – 500.0 mL, Graduation Interval: 5.0 mL, Pharmaceutical Grade, Double Scale

**DESCRIPTION:** This graduated cylinder meets ASTM E-1272 Class A, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

**TECHNICAL SPECIFICATIONS:**

Capacity (mL)	500
Approximate O.D. x Height (mm)	53 x 362
Graduation Range (mL)	25.0 – 500.0
Graduation Interval (mL)	5.0
Tolerance (mL)	± 2.0