



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9854



CYLINDER, GRADUATED GLASS TYPE I, CLASS B (250 mL)

Graduated Range: 10.0 – 250.0 mL, Graduation Interval: 2.0 mL, Pharmaceutical Grade, Double Scale

DESCRIPTION: This graduated cylinder meets ASTM E-1272 Class B, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	250
Approximate O.D. x Height (mm)	41 x 301
Graduation Range (mL)	10.0 – 250.0
Graduation Interval (mL)	2.0
Tolerance (mL)	± 1.40