



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9855



CYLINDER, GRADUATED GLASS TYPE I, CLASS A (100 mL)

Graduated Range: 1.4 – 25.0 mL, Graduation Interval: 0.2 mL, Pharmaceutical Grade, Double Scale

DESCRIPTION: This graduated cylinder meets ASTM E-1272 Class A, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	100
Approximate O.D. x Height (mm)	30 x 236
Graduation Range (mL)	5.0 – 100.0
Graduation Interval (mL)	1.0
Tolerance (mL)	± 0.40