



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9856



CYLINDER, GRADUATED GLASS TYPE I, CLASS A (25 mL)

Graduated Range: 1.4 – 25.0 mL, Graduation Interval: 0.2 mL, Pharmaceutical Grade, Single Scale

DESCRIPTION: This graduated cylinder meets ASTM E-1272 Class A, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

TECHNICAL SPECIFICATIONS:

Capacity (mL)	25
Approximate O.D. x Height (mm)	17 x 210
Graduation Range (mL)	1.4 – 25.0
Graduation Interval (mL)	0.2
Tolerance (mL)	0.15