



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9858**



**CYLINDER, GRADUATED GLASS TYPE I, CLASS A (5 mL)**

Graduated Range: 0.5 – 5 mL, Graduation Interval: 0.1 mL, Pharmaceutical Grade, Single Scale

**DESCRIPTION:** This graduated cylinder meets ASTM E-1272 Class A, Type I Specifications & USP standards for volumetric glassware. Hexagonal base provides extra stability and all markings are in durable, white enamel. Comes with a protective, polyethylene, hexagonal, bumper guard.

**TECHNICAL SPECIFICATIONS:**

Capacity (mL)	5
Approximate O.D. x Height (mm)	7 x 126
Graduation Range (mL)	0.5 – 5
Graduation Interval (mL)	0.1
Tolerance (mL)	0.05