



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9863



CONICAL CYLINDER, GRADUATED GLASS TYPE I (50 mL)

Graduated Range: 10 – 50 mL, Graduation Interval: 5 mL, Pharmaceutical Grade, Dual Scale

DESCRIPTION: This conical cylinder is used to perform all sort of measurements.

FEATURES:

- ▣ Scales in both metric and U.S. customary units
- ▣ "TD" appears on the cylinder and indicates it is calibrated to deliver
- ▣ Designed from ASTM Specification E1094, Type 2 requirements (also NBS Handbook 44)
- ▣ Manufactured from 33-expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

TECHNICAL SPECIFICATIONS:

Capacity (mL)	50
Height (mm)	139
Graduation Range (mL)	10 – 50
Graduation Interval (mL)	5
Tolerance (mL)	Up to 25 ± 0.6, above 25 ± 0.8