



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 9959**



### **ZIP-LOCK BAG (Amber, 3-mil thick)**

**6" × 8"**

**DESCRIPTION:** This zip-lock bag is made of translucent polyethylene material and it provides UV protection for light-sensitive items. The side-weld construction and zipper seal helps to protect from dirt, moisture, and other contaminants. The usable dimension is slightly smaller.

Meets requirements of the U.S. Pharmacopoeia, Revision 35 regarding protection against light transmission.

#### **TECHNICAL SPECIFICATIONS:**

|                       |              |
|-----------------------|--------------|
| Dimension, W × D (in) | 6 × 8        |
| Thickness (in)        | 0.003        |
| Color                 | Amber        |
| Material              | Polyethylene |