



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 9986



VAGINAL CREAM APPLICATOR (4 mL)

Plastic, Disposable, Natural

DESCRIPTION: This applicator can be easily loaded by attaching it directly to a #16 ointment tube. It has a unique feature which provides a stop at a specific location to accommodate the user's unique single dose requirement which prevents the user from administering an improper dose.

TECHNICAL SPECIFICATIONS:

Capacity	4 mL
Dosage Increment	1 g
Color	Natural